Macroplastique® Implants

*For the Treatment of Stress Urinary Incontinence
and Vesicoureteral Reflux*

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Lot Number

Manufacturer

Authorized Representative in a European Community

INDICATIONS

Stress Urinary Incontinence: Macroplastique Implants are indicated for the treatment of stress urinary incontinence (SUI) resulting from an incompetent sphincteric mechanism. Macroplastique may be used as a primary treatment or as a secondary treatment following previous operative procedures. When injected between the mid-urethra and bladder neck, this urethral bulking agent aids in the restoration of urinary continence.

Vesicoureteral Reflux: Macroplastique Implants are indicated for the treatment of vesicoureteral reflux (VUR). When injected at the ureteral orifice, the tissue bulking provides occlusive pressure to reduce or eliminate vesicoureteral reflux.

DESCRIPTION

Macroplastique is a soft-tissue bulking agent consisting of solid, irregularly textured, medical grade polydimethylsiloxane elastomer implants suspended in a hydrogel carrier of polyvinylpyrrolidone. The implants are held in place at the injection site when the hydrogel carrier is replaced by body fluids and host fibroblasts subsequently deposit collagen around the implants. The hydrogel is removed by the reticuloendothelial system and excreted unmetabolized from the body through the kidneys.

CONTRAINDICATIONS

- Acute urogenital tract inflammation
- Gross utero-vaginal prolapse
- Untreated detrusor instability
- Detrusor hyperreflexia
- Neuropathic bladder
- Overflow incontinence
- Pregnancy
- Within one year post-partum
- Within 12 weeks of a previous Macroplastique treatment
- Within 12 weeks of a previous mid-urethral sling placement

WARNINGS

- Improper injection technique and injected volume may result in unsatisfactory performance and require corrective surgery.
- Improper patient selection may result in unfavourable results.
- Collect and microbiologically evaluate urine to ensure the absence of urinary tract infection. Do not proceed if the patient has a urinary tract infection.

- Macroplastique is generally not recommended for placement in fibrotic tissue or in patients previously treated with polytetrafluoroethylene or rigid particulate implants for incontinence, or with malignancy.
- Injections into the urethra or the subureteral space require precision and should only be performed by surgeons trained in diagnostic and therapeutic cystoscopy.
- Do not use indwelling catheters post-operatively.
- If there is any visual damage to the syringe and/or the inner or outer package, the integrity of the product is not guaranteed and should not be used. Return the entire package to Uroplasty or your distributor for replacement.
- For safe and accurate injection of Macroplastique, only use needles supplied and/or approved by Uroplasty.
- Reuse of opened Macroplastique is associated with risks of infection and inability to express the product through a needle due to moisture loss.

ADVERSE REACTIONS

- General complications associated with anaesthesia, procedural methods, and patient's general tolerance to any injected material
- Dysuria, hematuria and frequency of micturition may occur postoperatively. Instruct the patient to contact the surgeon if any of these conditions persist past 48 hours.
- Post-operative retention may occur which may necessitate intermittent catheterization. If the patient remains unable to void freely, continued intermittent catheterization may be necessary.
- If the procedure is followed by reactions in excess of those expected, diagnostic and therapeutic procedures should be implemented at the discretion of the surgeon.
- Notify Uroplasty or your distributor of product-related complications occurring with the use of this device.

TECHNIQUE CONSIDERATIONS

- Surgeons should be familiar with the flow characteristics of Macroplastique prior to injection.
- The injected volume of Macroplastique is approximately equal to the permanent volume it will occupy when the collagen network is formed.
- Administer product slowly and wait 30 seconds before removing needle from injection site. If injection is performed too quickly, Macroplastique can continue to pass through the needle under pressure.
- Inject Macroplastique within the submucosal tissue of either the urethra (for SUI) or the subureteral space (for VUR). Use proper tunnelling technique to avoid injecting product too superficially or too deeply (see "Tunnelling Technique").
- Superficial injection will lead to mucosal damage and loss of the implants. Implantation below the submucosa will result in a bolus that may be too deep to produce the desired tissue bulking.
- Do not open the package until immediately prior to use; product may dry out and become difficult to inject.
- Follow strict aseptic technique during injection.

PRE-OPERATIVE PATIENT CARE

1. Collect and microbiologically evaluate urine to ensure absence of urinary tract infection. Rule out systemic infections.
2. Administer prophylactic broad-spectrum antibiotics.
3. Ask patient to void prior to procedure.
4. Place patient in a modified lithotomy position. Clean area with surgical skin preparation and position sterile drapes.
5. Prepare patient with local, conscious sedation or general anaesthesia.
6. When using local anaesthetic, insert 20 cc of 2% lidocaine gel into the urethra with the supplier's dispensing probe. Administer transurethral or periurethral lidocaine block. Wait for local anaesthetic to take effect. Note: individual anaesthetic injections >1 ml may cause hydro-dissection of urethral tissues resulting in an

unstable bolus deposition and may compromise the successful outcome of the procedure.

7. Fill bladder to 50% capacity.

PROCEDURE: SYSTEM PREPARATION

1. Aseptically remove the Macroplastique syringe from its packaging. Keep cap on tip of syringe.
2. Slide adapter over the syringe until it fits down around the syringe flanges. Lock the assembled adapter and syringe onto the Administration Device. Firmly grasp both the adapter and the rotating collar to secure the syringe to the Administration Device. The rotating collar holds the Administration Device in a constant position while the syringe/needle is rotated.
3. Aseptically remove the endoscopic needle from its packaging. Without removing the needle protective sleeve, attach the endoscopic needle by firmly tightening the needle hub onto the syringe Luer lock tip. Remove needle protective sleeve.
4. Prior to inserting the needle into the cystoscope, prime the needle by pulling and releasing the Administration Device lever several times until product is visible at the tip of the needle.
5. Press the release mechanism at the top of the Administration Device to stop product flow.

TREATMENT FOR MALE AND FEMALE INCONTINENCE

Note: When using the Macroplastique Implantation System (MIS) refer to its accompanying insert.

Macroplastique Injection Locations

Females	Inject 1.5 – 2.0 cm distal from the bladder neck
Males external sphincter functional length ≤ 1.5 cm	Inject at the most distal part of the external sphincter but before the bulbous urethra.
Males external sphincter functional length ≥ 1.5 cm	Inject 10 mm distal to the midpoint of the external sphincter. Area should have good vascularization and minimal scar tissue.

Approximate Treatment Volumes and Positions

Macroplastique should be injected in 2-3 positions around the urethra to achieve urethral coaptation. Injection positions and volumes should be adjusted according to the patient's history of incontinence surgery (i.e., bladder neck suspension, sling procedures, etc.) and morphology of the bladder neck and urethra.


	Approximate positions			Total volume
	6 o'clock	10 o'clock	2 o'clock	
Females - Primary treatment	≤ 2.5 ml	≤ 1.25 ml	≤ 1.25 ml	≤ 5 ml
Males - Primary treatment	≤ 2.5 ml	≤ 2.5 ml	≤ 2.5 ml	≤ 7.5 ml
Retreatment after previous urethral bulking procedure	Volume and position dependent on morphology of the bladder neck and urethra			≤ 5 ml
Secondary treatment after previous surgical procedure (i.e., sling)	Surgeon discretion	≤ 2.5 ml	≤ 2.5 ml	≤ 5 ml

Injection Procedure

1. Insert cystoscope into urethra and advance needle through working channel to visualize needle tip.
2. Advance needle with bevel facing the center of urethra.
3. To ensure appropriate placement of Macroplastique in the submucosal tissue, use the tissue tunneling technique during the procedure (Figures A-D). Wait approximately 30 seconds before withdrawing needle from tissue to limit product loss from injection site.
4. Locate the injection position within urethral lumen and tilt scope to a 30-45° angle (Figure A).

5. Insert needle tip into urethral tissue at this angle and then advance needle 0.5 cm in depth. Angle scope to 0°, parallel with the urethra (Figure B).
6. Advance needle another 0.5 cm to create a tissue tunnel (Figure C).
7. Inject a small amount of Macroplastique to confirm correct needle placement within mucosa (Figure D). If needle is properly placed, tissue bulking in the urethral mucosa should immediately be visible with the cystoscope. If it does not appear, withdraw and reposition needle more superficially and inject again.

Figures A-D: Tissue Tunneling Technique

8. Inject product slowly. Wait a few seconds between each pull of the Administration Device lever.
9. Once desired tissue-bulking is achieved, repeat injection procedure in the other location(s) to achieve urethral coaptation.
10. To change syringe during procedure, press the release mechanism at the top of the Administration Device, retract dosage rod and remove syringe and needle. Reload system (see "Procedure: System Preparation").
11. Avoid passing cystoscope over injection site to prevent product disruption.

Post-operative Patient Care

- Use a small intermittent catheter (8-12 Fr.) to drain the bladder when necessary.
- Prescribe prophylactic broad-spectrum antibiotics.
- Upon discharge, patient must be able to void freely. In cases of urinary retention, intermittent catheterization is advised until normal voiding resumes. Indwelling urethral catheterization is not recommended.
- Patients should avoid all forms of sexual intercourse, rigorous exercise and heavy lifting for 1 month following surgery.

ADDITIONAL TREATMENT

If further treatment is required, wait 12 weeks to allow the host tissue to fully mature. In general, Macroplastique should be injected distal to the initial injection sites. However, injection sites should be adjusted according to morphology of the bladder neck and urethra.

TREATMENT FOR VESICoureTERAL REFLUX (VUR)

1. Introduce needle 2-3 mm below the refluxing ureteral orifice at the 6 o'clock position. While remaining in the submucosal layer, advance needle approximately 4-5 mm until needle tip is directly beneath the ureteral orifice.
2. Inject a small amount of Macroplastique to verify correct needle placement within the mucosa. If the needle is properly placed, tissue bulking will immediately appear at the ureteral orifice.
3. Inject 0.2 to 1.0 ml of Macroplastique into the submucosal layer of the subureteral space until a volcano-like appearance is achieved, resulting in a crescent-shaped ureteral opening located on top of the augmented tissue.
4. Wait approximately 30 seconds before withdrawing needle from tissue to limit product loss from the injection site.

Post-operative Patient Care

Prescribe prophylactic broad-spectrum antibiotics.

Additional Treatment

If further treatment is required, wait 12 weeks to allow the host tissue to fully mature. In general, Macroplastique should be injected distal to the initial injection site. The goal is to create a more distinct formation with a slit-like orifice.

PATIENT INFORMATION AND INFORMED CONSENT

Uroplasty relies on the surgeon to advise the patient of all potential risks and complications associated with the proposed surgical procedure, including providing a comparison of the risks, benefits and complications of alternate procedures.

DISCLAIMER OF WARRANTIES

Uroplasty excludes all warranties and responsibilities for:

- Improper use of the product, and/or
- Failure to follow instructions in this insert and Macroplastique Implants accessories.

ORDERING INFORMATION

Pre-filled syringes of Macroplastique are supplied sterile and latex-free for single use in two volumes:

- MPQ-1.5 Approximately 1.5 ml
- MPQ-2.5 Approximately 2.5 ml

Macroplastique accessories include:

- Administration Device
- Macroplastique Implantation System
- Flexible Endoscopic Needles
- Rigid Endoscopic Needles

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Macroplastique and Uroplasty are registered trademarks of Uroplasty, Inc.
Macroplastique und Uroplasty sind eingetragene Marken von Uroplasty, Inc.
Macroplastique y Uroplasty son marcas registradas de Uroplasty, Inc.
Macroplastique et Uroplasty sont des marques déposées de Uroplasty, Inc.
Macroplastique e Uroplasty sono marchi registrati di Uroplasty, Inc.
Macroplastique en Uroplasty zijn gedeponeerde handelsmerken van Uroplasty, Inc.